


Particulate Filter CF1 (lot de 2) HEPA H14

References and Filter Type

HP	H12111021
----	-----------

Dimensions and Weight

Filter Dimensions	190 x 403 x 61 mm
Protruding Seal(s)	0
Filter Weight	Env. 5 Kg



Filter Composition

HEPA H14 Filter	Metal structure
	Non-woven fiber mesh

Features

HEPA H14 Filter	Particulate filtration for particles/powders or fumes
-----------------	---

Filtration Efficiency

Particulate Filtration	99.995% efficiency according to MPPS method (Most Penetrating Particle Size) of EN1822-1:1998 standard
------------------------	--

Warnings

Your protection is only guaranteed after approval by the Erlab® laboratory.

Contact your supplier and benefit from the **ESP® (Erlab Safety Program)** to validate the compatibility of the filter type with your applications.

- This filter must be stored in its original packaging and in a dry place.
- It is recommended to wear gloves, safety glasses, and a lab coat when handling this filter.