


## Particulate Filter GF4 ULPA U16

References and Filter Type	
U16	H11074061

Dimensions and Weight	
Filter Dimensions	745.5 x 394.5 x 89.5 mm
Protruding Seal(s)	0
Filter Weight	Env. Kg



Filter Composition	
ULPA U16 Filter	Metal structure
	Non-woven fiber mesh

Features	
ULPA U16 Filter	Particulate filtration for particles/powders or fumes

Filtration Efficiency	
Particulate Filtration	99,99995% efficiency according to MPPS method (Most Penetrating Particle Size) of EN1822-1:1998 standard

Warnings	
<p><b>Your protection is only guaranteed after approval by the Erlab® laboratory.</b>            Contact your supplier and benefit from the <b>ESP® (Erlab Safety Program)</b> to validate the compatibility of the filter type with your applications.</p> <ul style="list-style-type: none"> <li>- This filter must be stored in its original packaging and in a dry place.</li> <li>- It is recommended to wear gloves, safety glasses, and a lab coat when handling this filter.</li> </ul>	