


## Particulate Filter F3 HEPA H14

### References and Filter Type

HP	H11013021
----	-----------

### Dimensions and Weight

Filter Dimensions	567.5 x 404 x 56 mm
Protruding Seal(s)	0
Filter Weight	Env. 4 Kg



### Filter Composition

HEPA H14 Filter	Metal structure
	Non-woven fiber mesh

### Features

HEPA H14 Filter	Particulate filtration for particles/powders or fumes
-----------------	---

### Filtration Efficiency

Particulate Filtration	99.995% efficiency according to MPPS method (Most Penetrating Particle Size) of EN1822-1:1998 standard
------------------------	--

### Warnings

**Your protection is only guaranteed after approval by the Erlab® laboratory.**

Contact your supplier and benefit from the **ESP® (Erlab Safety Program)** to validate the compatibility of the filter type with your applications.

- This filter must be stored in its original packaging and in a dry place.
- It is recommended to wear gloves, safety glasses, and a lab coat when handling this filter.